

## Troubleshooting .....

If problem cannot be corrected, consult your Dealer/Home Care equipment provider.

PROBLEM	PROBABLE CAUSE	REMEDY
A. No flow	<ol style="list-style-type: none"><li>1. Cylinder valve closed</li><li>2. Converter in "OFF" position</li><li>3. Cylinder empty</li><li>4. Dial not at setting</li><li>5. Converter not sensing breath</li></ol>	<ol style="list-style-type: none"><li>1. Turn on cylinder</li><li>2. Set to prescribed setting</li><li>3. Replace cylinder</li><li>4. Set dial so indicator points to a setting</li><li>5. Check position of cannula in nose Do not breathe through mouth</li></ol>
B. Leaking at Cylinder Connection	<ol style="list-style-type: none"><li>1. Missing or defective O-ring Seal</li></ol>	<ol style="list-style-type: none"><li>1. Replace O-ring Seal</li></ol>

**Precision Medical, Inc.**  
300 Held Drive • Northampton, PA 18067  
800-272-7285 • 610-262-6090 • fax: 610-262-6080

504442 (E) 4/05 Printed in USA


## Setting up your

# ***EASYPULSE5*** *Oxygen Conserving Regulator*

## **1900 Series (CGA 540 Connection)**




**PRECISION**  **MEDICAL**®

 Quick Start Guide only  
Refer to User Manual for complete operating instructions, warnings and cautions.

## EasyPulse5 Oxygen Conserving Regulator (Conserver)


# PARTS



## STEP 1

Inspect Conserver and Cylinder Valve for damage or contaminants. Do NOT use if any are found. “Crack” cylinder valve by opening the valve slightly and closing.

Cylinder Valve



## STEP 2

Ensure Conserver is in the OFF position and O-ring Seal is in place. Align Nipple with recessed hole in the cylinder. Hand tighten Knob. Do NOT use tools.

View of properly attached conserver



## STEP 3

When using a Precision Medical Carry Bag, slide the cylinder into the bag with the gauge facing the mesh. Attach the cannula to the Conserver’s Outlet Connection.



## STEP 4

Slowly open the cylinder valve. Turn Dial to the prescribed setting. Zip the bag closed.



## STEP 5

Place cannula in nostrils and breathe through nose to begin Oxygen flow. If using a bag adjust the straps for maximum comfort.



Cannula prongs in nostrils



Typical bag carrying configurations

## STEP 6

When finished using: Turn the cylinder valve to the OFF position. Turn Conserver Dial to “2 LPM Continuous” position. Wait for Oxygen to stop flowing from Conserver. Remove Conserver from Cylinder Valve. Store Conserver and Carry Bag in a safe place.