

PARTS

EasyPulse5 Oxygen Conserving Regulator (Conservor)

STEP 1 Inspect Conservor and Cylinder Valve for damage or contaminants. Do NOT use if any are found. "Crack" cylinder valve by opening the valve slightly and closing.

STEP 2 Ensure Conservor is in the OFF position and O-ring Seal is in place. Align Nipple with recessed hole in the cylinder. Hand tighten Knob. Do NOT use tools.

Setting up your EASYPULSE⁵ Oxygen Conserving Regulator Model 195405 (CGA 540 Connection)

PRECISION MEDICAL

Quick Start Guide only
Refer to User Manual for complete operating instructions, warnings and cautions.

TROUBLESHOOTING

If problem cannot be corrected, consult your Dealer/Home Care equipment provider.

PROBLEM	PROBABLE CAUSE	REMEDY
A. No flow	1. Cylinder valve closed 2. Conservor in "OFF" 3. Cylinder empty 4. Dial not at setting 5. Conservor not sensing breath	1. Turn on cylinder position setting 2. Set to prescribed position setting 3. Replace cylinder 4. Set dial so indicator points to a setting 5. Check position of cannula in nose Do not breathe through mouth 1. Replace O-ring Seal
B. Leaking at Cylinder Connection	1. Missing or defective O-ring Seal	1. Replace O-ring Seal

STEP 3 When using a Precision Medical Carry Bag, slide the cylinder into the bag with the gauge facing the mesh. Attach the cannula to the Conservor's Outlet Connection.

STEP 4 Slowly open the cylinder valve. Turn Dial to the prescribed setting. Zip the bag closed.

STEP 5 Place cannula in nostrils and breathe through nose to begin Oxygen flow. If using a bag adjust the straps for maximum comfort.

Cannula prongs in nostrils Typical bag carrying configurations

STEP 6 When finished using: Turn the cylinder valve to the OFF position. Turn Conservor Dial to "2 LPM Continuous" position. Wait for Oxygen to stop flowing from Conservor. Remove Conservor from Cylinder Valve. Store Conservor and Carry Bag in a safe place.