

# CHEMETRON STYLE QUICK CONNECTS

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 <b>1/8 NPT Male</b>	3101	3201	3301	3401	3601
 <b>1/8 NPT Female</b>	3102	3202	3302	3402	3602
 <b>1/4 NPT Male</b>	3103	3203	3303	3403	3603
 <b>DISS Male</b>	3105	3205	3305	3405	3605
 <b>DISS Hex Nut</b>	3106	3206	3306	3406	3606
 <b>DISS Hand Tight</b>	3108	3208	3308	3408	3608
 <b>Hose Barb</b>	<b>1/4 Hose Barb</b> 3111	<b>1/4 Hose Barb</b> 3211	<b>1/4 Hose Barb</b> 3311 <b>5/16 Hose Barb</b> 3312	<b>1/4 Hose Barb</b> 3411	<b>1/4 Hose Barb</b> 3611 <b>5/16 Hose Barb</b> 3615
 <b>DISS Male with Check Valve</b>	3114	3214	3314	3414	3614