












MEDSTAR STYLE QUICK CONNECTS

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum
 1/8 NPT Male	8101	8201	8301
 1/8 NPT Female	8102	8202	8302
 1/4 NPT Male	8103	8203	8303
 DISS Male	8105	8205	8305
 DISS Hex Nut	8106	8206	8306
 DISS Hand Tight	8108	8208	8308
 Hose Barb	1/4 Hose Barb 8111	1/4 Hose Barb 8211	1/4 Hose Barb 8311
 DISS Male with Check Valve	8114MS	8214	8314