

PURITAN-BENNETT STYLE QUICK CONNECTS

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Male	4101	4201	4301	4401	4601
 1/8 NPT Female	4102	4202	4302	4402	4602
 1/4 NPT Male	4103	4203	4303	4403	4603
 DISS Male with Knob	4105	4205	4305	4405	4605
 DISS Hex Nut	4106	4206	4306	4406	4606
 DISS Hand Tight	4108	4208	4308	4408	4608
 Hose Barb with Knob	1/4 Hose Barb 4111	1/4 Hose Barb 4211	1/4 Hose Barb 4311 5/16 Hose Barb 4312	1/4 Hose Barb 4411	1/4 Hose Barb 4611 5/16 Hose Barb 4615
 DISS Male with Knob and Check Valve	4114	4214	4314	4414	4614