













SCHRADER STYLE QUICK CONNECTS

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 Nitrogen
 1/8 NPT Male	6101	6201	6301	6401	6501
 1/8 NPT Female	6102	6202	6302	6402	6502
 1/4 NPT Male	6103	6203	6303	6403	6503
 DISS Male with Knob	6105	6205	6305	6405	6504
 DISS Hex Nut	6106	6206	6306	6406	6506
 DISS Hand Tight	6108	6208	6308	6408	6508
 Hose Barb	1/4 Hose Barb 6111	1/4 Hose Barb 6211	1/4 Hose Barb 6311	1/4 Hose Barb 6411	1/4 Hose Barb 6511
 DISS Male with Knob and Check Valve	6114	6214	6314	6414	6514

Schrader Piston Identification	
	Hex Piston (standard)
	Swivel Piston (add "S" after part number)