

PEDIATRIC FLOWMETERS



Extreme accuracy when it counts.

Glass Flowtube....

The glass flowtube assures consistent readings and allows superb accuracy for precise low flow requirements.

“Stork” Safety Label....

With the glance of an eye, the stork label quickly identifies the low flow flowmeter, preventing accidental use of wrong flowmeter in neonatal application.

Shock Protection....

The spring above the flowtube and the rubber flow tube bushing absorbs shock to protect the glass tube from breakage.

Warranty....

All flowmeters are backed with a five (5) year warranty on the needle valve and a lifetime warranty on the housing.

Black Glass Float Indicator....

The black glass float ball enhances visibility in order to assure accurate flow settings.

Durable Housing....

The polycarbonate housing is strong and durable. It is impervious to chemical solutions normally used in a clinical setting.

Front Flow Control Knob....

The flow control knob is on the front of the flowmeter in clear sight of the operator. Side mount knobs are susceptible to damage during installation of additional equipment.

ORDERING INFORMATION (PEDIATRIC FLOWMETERS)



O₂
0-200 CC



O₂
0-1 LPM



O₂
0-3 LPM

	No Adaptor	6MFA1001	4MFA1001	1MFA3001
	DISS Female Hex Nut	6MFA1002	4MFA1002	1MFA3002
	DISS Female Hand Tight	6MFA1003	4MFA1003	1MFA3003
	DISS Male	6MFA1004	4MFA1004	1MFA3004
	Ohmeda	6MFA1005	4MFA1005	1MFA3005
	Chemetron	6MFA1006	4MFA1006	1MFA3006
	Oxequip	6MFA1007	4MFA1007	1MFA3007
	Puritan-Bennett	6MFA1008	4MFA1008	1MFA3008
	Schrader	6MFA1009	4MFA1009	1MFA3009

SPECIFICATIONS

Model	6MFA1001	4MFA1001	1MFA3001
Flow Range	0-200 cc	0-1 lpm	0-3 lpm
Gas	Oxygen	Oxygen	Oxygen
Increments	20 cc starts at 20 cc	.1 lpm starts at .1 lpm	.125 lpm from .125 to 1 lpm .25 lpm from 1 to 3 lpm
Accuracy	± 10 cc from 0-100 cc ± 14 cc from 100-200 cc	± .05 lpm	± .15 lpm
Max Flush Flow	500 cc	5 lpm	40 lpm

Flowmeters calibrated at 50 psi (.34 bar), 70°F (21°C), standard atmospheric pressure.
Specifications are subject to change without prior notice.



ISO 9001
Certified

Precision Medical, Inc.
300 Held Drive • Northampton, PA 18067
Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285
Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240