

COMPACT FLOWMETERS



State of the art design!

They're Lightweight & Durable...

The polycarbonate housing extends over the lightweight aluminum body providing superior durability. The housing is impervious to chemical solutions normally used in a clinical setting.

They're Easy to Read...

The state of the art housing incorporates a white background that makes the flow settings easier to read.

They're Compact...

The contemporary design results in a smaller unit that fits more easily into space-cramped wall rails. Unlike side mounted knobs which are susceptible to damage, the flow control knobs are on the front of the flowmeter in clear sight of the operator.

They're Accurate...

Our flowmeters use black glass flow indicator balls. The properties of black glass, superior sphericity and exceptionally uniform density, insure that the flowmeter is more accurate than flowmeters utilizing plastic floats.

They're Dependable...

The flow tube has a lifetime warranty and the needle valve carries a five (5) year warranty.



SPECIFICATIONS

Model	Gas	Scale	Increments	Accuracy
0-8 lpm	Oxygen	0-8 lpm	.5 lpm from .5 to 8 lpm	±.25 lpm
0-15 lpm	Oxygen	0-15 lpm	.5 lpm from .5 to 5 lpm 1 lpm from 5 to 15 lpm	±.25 lpm from .5 to 5 lpm ±.5 lpm above 6 to 15 lpm
0-15 lpm	Air	0-15 lpm	.5 lpm from .5 to 5 lpm 1 lpm from 5 to 15 lpm	±.25 lpm from .5 to 5 lpm ±.5 lpm above 6 to 15 lpm

Flowmeters calibrated at 50 PSI, 70°F, standard atmospheric pressure. Specifications are subject to change without prior notice.

ORDERING INFORMATION

		O2 0-8 LPM	O2 0-15 LPM	Air 0-15 LPM
	No Adaptor (1/8 NPTF)	8MFA4001	8MFA1001	8MFA2001
	DISS Female Hex Nut	8MFA4002	8MFA1002	8MFA2002
	DISS Female Hand Tight	8MFA4003	8MFA1003	8MFA2003
	DISS Male	8MFA4004	8MFA1004	8MFA2004
	Ohmeda	8MFA4005	8MFA1005	8MFA2005
	Chemetron	8MFA4006	8MFA1006	8MFA2006
	Oxequip	8MFA4007	8MFA1007	8MFA2007
	Puritan-Bennett	8MFA4008	8MFA1008	8MFA2008
	Schrader	8MFA4009	8MFA1009	8MFA2009



ISO 9001
Certified

Precision Medical, Inc.
300 Held Drive • Northampton, PA 18067
Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285
Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240
Visit us on the Web: www.precisionmedical.com

503656rev2 6/05 7.5M Printed in USA